

The background features a series of overlapping, curved lines in shades of blue and red, creating a sense of motion and depth. The lines are thick and have a slight gradient, giving them a three-dimensional appearance. They are arranged in a roughly circular pattern, with some lines crossing each other.

# SEEBECK MEASUREMENT SYSTEM

## Thermoelectric Measurements

# THE SEEBECK MEASUREMENT SYSTEM




# THE SEEBECK EFFECT

- The Seebeck Effect, or thermoelectric effect
  - Direct conversion of temperature differences into electricity
- Discovered in 1821 by Thomas Johann Seebeck
  - Found that a voltage existed between two ends of a metal bar when a temperature gradient,  $\Delta T$ , existed in the bar
  - When a closed loop is made of two metals with a temperature difference, a magnetic field is produced.
    - Effect is that a voltage, or thermoelectric EMF, is created in the presence of two metals with a temperature difference between the junctions


# THE SEEBECK EFFECT

$$dV = S_{AB}(K_h - K_c)$$



- A and B are a pair of dissimilar metals, whose two junctions are held at different temperatures
- $K_h - K_c$  is the difference in temperature of the hot and cold junctions
- $S_{AB}$  is the relative Seebeck Coefficient, thermoelectric power
  - Varies with the level of the temperature at which the temperature difference occurs
- $dV$  – voltage difference across the terminals of an open circuit
  - Does not depend on the distribution of temperature along the metals between the junctions

# THE SEEBECK EFFECT



$$V = \int_{T_1}^{T_2} (S_B(T) - S_A(T)) dT$$

$$V = (S_B - S_A)(T_2 - T_1)$$

Simplify the equation if the Seebeck Coefficients are effectively constant of the measured temperature range

- If the circuit is closed, a current will flow through in the metals, detected by:
  - Magnetic field produced around the wires
  - Joule heating produced by the resistance in the wires
  - Galvanometer or ammeter placed in circuit to measure the current

# CAUSES OF THE SEEBECK EFFECT



## ➤ Charge Carrier Diffusion

- Charge carriers diffuse when one end of a conductor is at a different temperature than the other
  - Hot carriers diffuse from the hot end to the cold end
  - Cold carriers diffuse from the cold end to the hot end
- Motion of charge carriers results in an electrical current


# CAUSES OF THE SEEBECK EFFECT

## ➤ Phonon Drag

- A phonon is a quantized mode of vibration occurring in a rigid crystal lattice
  - Not always in local thermal equilibrium – they move along the thermal gradient
- Loss of momentum by interaction with carriers and imperfections in the lattice
  - If the phonon-electron interaction is predominant, phonons tend to push electrons to one end of the material
    - Occurs strongly at the Debye temperature
    - Region of thermopower vs. temperature is highly variable under a magnetic field

# TAKING MEASUREMENTS

- Two pairs of thermocouples:
  - One based on junctions of Cu-Metal with known properties
  - One based on Cu-Metal of unknown properties
- Computer controlled heater, located close to working junctions, far from reference




# TAKING MEASUREMENTS

- Principle of operation:
  - Assume all four thermocouples have same temperature
    - $V_1$  and  $V_2 = \text{zero}$
  - Apply power to heater to create temperature difference between working and reference junctions
    - Get some non-zero values for  $V_1$  and  $V_2$
    - Assume the temperature difference is the same for both pairs because the sample stage is symmetrical
      - Gives the ratio of the specific thermovoltages equal to the differences in the voltages  $V_1$  and  $V_2$

# TAKING MEASUREMENTS

- Difference is very small, so direct measurement will not give high accuracy
  - Eliminate the inaccuracy by taking measurements at two temperature points, and operate with the difference signal
  - Additionally, take multiple measurements at each point and average the results



Reference Wire

Sample "Wire"

# COMPONENTS IN A SEEBECK SYSTEM

- Pure high-pressure gas (greater than 1800 psi)
  - Gas Lines, Pressure Gauge, etc
- Filter/Dryer Apparatus
- Refrigerator
- Computer
- Temperature controller
- Circuit Breakout Box
- Vacuum Pump
- Seebeck Vacuum Chamber
- Seebeck Electronics and Software
- Sample Mounting Stage

# GAS, LINES, GAUGES, ETC

---


- 99.998% Pre-purified Nitrogen at 2640 psi or higher
- High Pressure Nitrogen Regulator
- High Pressure Nitrogen Lines (supplied)
- Gas Flow Meter (supplied)

# FILTER/DRYER SYSTEMS



# THE JOULE-THOMSON REFRIGERATORS

---




R2500-XX

# COMPUTER SYSTEM

---

- Pentium Processor, 1 GHz minimum
- CD-ROM Drive (4x)
- 250 MB RAM
- 50 MB free on hard drive
- 2 RS232 Serial Ports or USB Ports

# VACUUM PUMP AND ACCESSORY KIT



# SEEBECK ELECTRONICS





# SEEBECK VACUUM CHAMBER



# SEEBECK SAMPLE MOUNTING STAGE

For Temperatures  
below 400 K



For Temperatures  
above 400 K

# POSSIBLE TEMPERATURE RANGES

## Kelvin Scale

- 70 K to 580 K
- 80 K to 580 K
- 70 K to 730 K
- 80 K to 730 K
- 300 K to 730 K

## Centigrade Scale

- - 200 °C to 305 °C
- - 190 °C to 305 °C
- - 200 °C to 455 °C
- - 190 °C to 455 °C
- 25 °C to 455 °C

$$K = ^\circ C + 273$$

# WHY CARE ABOUT THE SEEBECK EFFECT?

---

- Automobiles and Fuel Consumption
- Energy Efficiency
- Alternate Energy Sources
- Thermoelectric refrigeration

# POTENTIAL END-USERS FOR SEEBECK

---

- Semiconductor industry
- Universities
- Energy and environmental engineers
- Thermoelectric Measurement
- Refrigeration industry